

Sequence File Syntax

- File extension - **.lseq**
- Preferred text editor – Sublime Text

General Format

- Lines starting with **#** are considered as comments. No action will be taken on these.
- Leaving blank lines is allowed.
- Avoid any tabs or indentation in the scripts.
- All keywords must be in FULL CAPS.

- All actionable commands have the following syntax –
`<action> <variable> <param1, param2, ...>`

- Only spaces must be used to separate words (and numbers) in actionable commands. No other punctuation should be used.

Commands Description

- Each **action** accepts certain **variables**. *Params* list depends on the **action** and **variable** chosen.
- Each **action** is described in the following subsections.

SET

- Drives a **variable** to a setpoint.
- Variables – **TEMP, FIELD, THETA, CURRENT**
- Params - `<setpoint> RATE <rate>`
- Examples –
 - `SET TEMP 300 RATE 6`
 - Set temperature to **300 K** at a rate of **6 K/min**
 - Do not exceed a rate of 10 K/min
 - `SET FIELD 1000 RATE 100`
 - Set field to **1000 Oe** at a rate of **100 Oe/s**
 - Never exceed a rate of 200 Oe/s
 - `SET THETA 90 RATE 1`
 - Set theta to **90 degrees** at a rate of **1 degree/s**
 - Never exceed a rate of 1 degree/s
 - `SET CURRENT 1.5e-6`
 - Set the current to **1.5×10⁻⁶ A**
 - Note that decimal as well as exponent representation is allowed for current setpoint.

WAIT

- Halts the sequence execution until setpoint of the variable is reached. Additionally delays proceeding of the sequence by specified time.
- Variables - **TEMP, FIELD, THETA, MEASURE**
- Params - *<delay time>*
- Examples –
 - **WAIT TEMP 600**
 - Waits until the temperature setpoint is reached. Then waits additional 600 **seconds** before proceeding to next lines.
 - Similar operation for **FIELD** and **THETA**
 - **WAIT MEASURE 600**
 - There is not setpoint to be reached here. This is just a provision to ensure that the sequence halts for a while until all instruments are properly initialized and ready to start meaningful measurements.
 - Halting at least 20 seconds after starting any measurement is recommended.

SCAN – END SCAN

- Starts a loop (something like a for loop) over a variable.
- Variables - **TEMP, FIELD, THETA, CURRENT**
- Params - *<start point> <stop point> <number of points> RATE <rate>*
- Examples –
 - **SCAN FIELD 0 90000 10 RATE 100**
WAIT FIELD 10
Additional steps to be executed in the loop
END SCAN
 - Start a scan over field. Fields setpoints loop over 0, 10000 Oe, ..., 90000 Oe (10 points, linearly spaced, from 0 to 90000).
 - **WAIT** command right after the **SCAN** command is recommended, to allow reaching the setpoint of the current loop.
 - Do not forget the **END SCAN** command at the end of the loop.
- Note – Nested loops are NOT supported.

START

- Start a measurement.
- Variables – **MEASURE, IV**
- Params – Not applicable
- Examples –
 - **START MEASURE**

- Starts continuous measurements based on the hardware settings (AC/DC system). Sequence execution proceeds to the next lines immediately. Measurements continue right until a STOP MEASURE command appears.
- `START IV`
 - Start an IV measurement, based on settings specified in the program. The IV plotting sub-VI pops up and closes after the IV measurement is complete. Sequence proceeds only after the IV measurement is complete.

STOP

- Stop the measurements started with the START MEASURE command.
- Variables – `MEASURE`
- Params – Not applicable
- Example –
 - `STOP MEASURE`

NEWFILE

- Create a new file for data accumulation. All data is appended in this file, over multiple START MEASURE – STOP MEASURE instances, until another file is created. File created has the chosen extension (default is .csv). A timestamp is also appended in the filename before the extension, in the format `_MMDD_hhmm`.
- Variables – Not applicable
- Params – *<filename without extension>*
- Examples –
 - `NEWFILE test_datafile`
 - Creates a new datafile with the name `test_datafile_1231_2359.csv` assuming .csv extension is chosen, and file was created on December 31 at 11:59 PM.